




EXPERIMENT-1

AIM — To estimate the amount of iron present in a given solution of Ferosus Ammonium Sulphate using potassium permanganate solution as intermediate.

MATERIAL REQUIRED — Mohr's Salt, KMnO_4 soln, dil. H_2SO_4 , Bullette, Funnel, Stand, White tile, Pipette, Volumetric Flask, Conical Flask, Heating Flask.

THEORY — Mohr's salt act as reducing agent and KMnO_4 act as O.A; so, the reaction between Mohr's Salt and KMnO_4 is a redox reaction. In this redox reaction ferrous ions from Mohr salt get oxidized and pink colour of Manganese present KMnO_4 which is the O_7 .

Oxidation state gets reduced to colourless Mn^{2+} state. KMnO_4 is powerful O.A. It helps in estimation of reducing agent like Mohr salt in acidic medium according to —



This titration is based upon oxidation-reduction titration. When Mohr's salt is

titrated against KMnO_4 in the presence of acidic medium by H_2SO_4 .

- PROCEDURE - 1) 10 ml of known Mohr's salt (11.92 g) was taken in a conical flask.
- 2) About half the test tube of H_2SO_4 was added to it.
- 3) The burette was filled with KMnO_4 soln.
- 4) The conical flask contents were titrated against.
- 5) 3 concordant reading were recorded.
- 6) The same process was repeated for an unknown soln. of Mohr's salt.

RESULT - The molarity of given Mohr's salt is 0.03 moles/l and strength is 11.92 g/l.

- PRECAUTIONS - 1) KMnO_4 is dark in colour so always read upper meniscus.
- 2) Rinse the pipette and burette before using.
- 3) Clean the apparatus with distilled water.
- 4) Strength of unknown soln should be taken upto 2 decimal places.

A. Observation table for Known Soln :

S No.	Vol. of Mohr's Salt soln (ml)	Burette Reading (ml)			Concordant Reading	Indicator
		Initial	Final	Volm of Soln		
1.	10	0	9.1	9.1	{ 9.1	KMnO ₄
2.	10	10	19.1	9.1	{ 9.1	(self-ind)
3.	10	20	29.1	9.1	{ 9.1	

B. Observation table for unknown soln :

S No.	Volm of Mohr's salt	Burette Reading (ml)			Concordant Reading	Indicator
		Initial	Final	Volm Used		
1.	10	0	12.5	12.5	{ 12.5	KMnO ₄
2.	10	12.5	25	12.5	{ 12.5	(self-ind)
3.	10	25	37.5	12.5	{ 12.5	

CALCULATIONS -

A. Normality of unknown Mohr's salt -

$$N_1 V_1 = N_2 V_2$$

$$\Rightarrow \frac{1}{40} \times 10 = N_2 \times 90$$

$$\Rightarrow N_2 = \frac{1}{36.4} N_1$$



$$N_2 V_2 = N_3 V_3 \Rightarrow \frac{1}{36.4} \times 12.5 = N_3 \times 10$$

EXPERIMENT - 2



AIM - To determine the strength of Potassium Dichromate soln using Ferric Ammonium Sulphate soln. as an intermediate.

MATERIAL REQUIRED - Mohr Salt, n-Phenyl Anthranilic Acid, Potassium Dichromate, dil H_2SO_4 , Bullette, Funnel, Stand, Measuring Flask, Pipette, Weighing bottle, Conical Flask, Volumetric Flask

THEORY - Acidic $K_2Cr_2O_7$ is a strong O.A., when it is added to Ferric Ammonium Sulphate soln., containing dil. H_2SO_4 only $FeSO_4$ is oxidised. $(NH_4)_2SO_4$ remain unchanged hence not shown in the eqn. The reaction takes place as follows -



complete reaction is -



- PROCEDURE - 1) 10 ml of Mohr's salt was taken in conical flask.
- 2) About half test tube of dil. H_2SO_4 was added to the conical flask.
- 3) A few drops of n-phenyl atraumatic acid was added as indicator.
- 4) The burette was filled with potassium dichromate.
- 5) The conical flask content was titrated against known $K_2Cr_2O_7$ soln. taken in burette.
- 6) The colour first turned green. On further adding $K_2Cr_2O_7$ from the burette, the colour of the soln. turned greenish purple. This was taken as end pt.
- 7) 3 concordant reading were taken.

RESULT - The strength of unknown $K_2Cr_2O_7$ soln. is ~~was given~~ 2.45 g/L.

- PRECAUTION - 1) Rinse the pipette and burette before use.
- 2) Don't use rubber cork.
- 3) Take accurate reading and don't go with avg. reading.

A. Observation Table for unknown $K_2Cr_2O_7$ soln -

Sl. No.	Vol. of Mohr Salt	Burette Reading (ml)			Concordanat Reading	Indicator used
		Initial	Final	Volm. used		
1.	10	10	16.5	6.5	6.5	n-Phenyl
2.	10	16.5	33	16.5	16.5	Anthracitic
3.	10	18.3	49.5	31.2	31.2	Acid.

B. Observation Table for unknown $K_2Cr_2O_7$ soln:

Sl. No.	Volm of	Burette Reading (ml)			Concordanat Reading	Indicator used
		Mohr Salt	Initial	Final	Volm used	
1.	10	9	18	8	8	n-Phenyl
2.	10	8	16	8	8	Anthracitic
3.	10	16	24	8	8	Acid.

CALCULATIONS —

$$N_1 V_1 = N_2 V_2 \Rightarrow \frac{1}{40} \times 16.5 = N_2 \times 10$$

$$\Rightarrow N_2 = 0.041 N$$

$$N_2 V_2 = N_3 V_3 \Rightarrow 0.04 \times 10 = N_3 \times 8$$

$$\Rightarrow N_3 = 0.05 N$$

∴ Normality of unknown $K_2Cr_2O_7$ soln is 0.05 N

→ Strength of unknown $K_2Cr_2O_7$ soln —

$$\text{Strength} = \text{Normality} \times \text{Eq. Mass}$$

$$= 0.05 \times \frac{291}{6}$$



$$= 2.45 \text{ g/L}$$

EXPERIMENT - 3

AIM - To determine the strength of a given CuSO_4 soln iodometrically.

MATERIALS REQUIRED - CuSO_4 soln, KI, Sodium Thiosulphate, Burette, Measuring Flack, Stand, Volumetric Flack, conical Flack, White Tile, Test tubes and funnel.

THEORY - When KI is added to a soln of CuSO_4 , a white cuprous iodide (Cu_2I_2) is precipitated and an equivalent amount of iodide is liberated. The free iodide is titrated with standard soln of hypo, using starch as an indicator. As soon as all the liberated iodide has been reduced to NaI , the blue colour of Iodo-starch complex, will disappear and the colour of ppt. in conical flack will be white.



PROCEDURE - 1) 10 ml of known CuSO_4 soln. was taken in a conical flack.
2) About half test tube of KI soln. was added to it.

- 3) The burette was filled with Sodium Thiosulphate soln.
- 4) Conical flask contents were titrated against the thiosulphate soln. taken in the
- 5) At this point, 5-6 drop of starch indicator was added to it.
- 6) The titration was continued with Sodium Thiosulphate soln, until the blue colour of complex just disappeared.
- 7) This was taken as the end point. Three concordant reading were taken.
- 8) The process was repeated for unknown CuSO_4 soln.

RESULT — The strength of CuSO_4 soln is
9.32 g/L.

PRECAUTION — 1) Iodine soln. attack rubber,
∴ it should always be taken in glass
stoppered burette.

2) Clean the apparatus with distilled water.

A. Observation Table for Known CUSO₄ Soln.

S. No.	Voln of CUSO ₄ (ml)	Burette Reading (ml)			Concordan Reading	Indicator Used
		Initial	Final	Total		
1.	10	0	14.2	14.2	14.2	Starch
2.	10	15	29.2	14.2	14.2	Starch
3.	10	30	44.2	14.2	14.2	Starch

RESULTS - 1. 10 ml CUSO₄ + 10 ml starch - 14.2 ml
 2. 10 ml CUSO₄ + 10 ml starch - 29.2 ml starch
 3. 10 ml CUSO₄ + 10 ml starch - 44.2 ml starch

S. No.	Voln of CUSO ₄ (ml)	Burette Reading (ml)			Concordan Reading	Indicator Used
		Initial	Final	Total		
1.	10	0	9.5	9.5	9.5	Starch
2.	10	10	9.5	9.5	9.5	Starch
3.	10	20	29.5	9.5	9.5	Starch

RESULTS - 1. 10 ml CUSO₄ + 10 ml starch - 9.5 ml starch
 2. 10 ml CUSO₄ + 10 ml starch - 9.5 ml starch
 3. 10 ml CUSO₄ + 10 ml starch - 29.5 ml starch

CALCULATIONS -

$$N_1 V_1 = N_2 V_2 \Rightarrow \frac{1}{N_2} \times 10 = N_2 \times 9.5 \Rightarrow N_2 = 0.025 \text{ N}$$

$$\Rightarrow N_2 = 0.025 \text{ N}$$

$$N_2 V_2 = N_3 V_3 \Rightarrow N_3 \times 10 = 0.025 \times 29.5 \Rightarrow N_3 = 0.027 \text{ N}$$

Take 1st 2nd soln. + 0.025 N

$$\text{Strength} = N \times \text{wt.}$$

$$= 0.027 \times 249.5$$

$$\text{Strength} = 9.32 \text{ g/l}$$

Strength = 9.32 g/l

Strength = 9.32 g/l

EXPERIMENT - 4

AIM - To calculate the concentration of unknown KMnO₄ solution using UV visible spectrometer.

MATERIALS REQUIRED - Test tubes, KMnO₄ soln of different concentrations, Light source, Monochromator, Sample, detector

THEORY - • Beer Lambert's Law - (Also called Beer Law) It is a relationship between the attenuation of light through a substance and the properties of that substance.

• Consider monochromatic light transmitted through a solution, with an incident intensity of I_0 and a transmitted intensity of I .

Then the transmittance $T = \frac{I}{I_0}$

$$T(\%) = 100 \times \frac{I}{I_0}$$

$$\text{The absorbance, } A = \log_{10} \frac{I_0}{I} = -\log_{10} T$$

• The Beer-Lambert law is linear relationship b/w the absorbance and the concentration, molar absorption coefficient and optical length.

$$A = \epsilon cl$$

$A \equiv$ Absorbance

$\epsilon \equiv$ Molar absorption coeff. ($M^{-1} cm^{-1}$)

$c \equiv$ Molar conc. (M); $l \equiv$ optical length

- PROCEDURE - 1) Prepare KMnO₄ stock solution ($\lambda_{max} = 525$ nm)
- 2) Prepare five standard soln of KMnO₄.
 - 3) Measure their absorbances at 525 nm.
 - 4) Plot calibration curve using the 5 absorbance and conc. values.
 - 5) calculate 'ε' [slope = $\epsilon \times d$]
 - 6) Measure absorbance of unknown sample.
 - 7) Using $A = \epsilon \cdot d$, calculate unknown KMnO₄ concentration.

RESULT - The concentration of unknown KMnO₄ soln is 0.344 M.

- PRECAUTIONS - 1) Clean the apparatus using distilled water.
- 2) Measure the readings carefully and with precision.

CALCULATION - Unknown Absorbance = 0.64


$$A = \epsilon cl \quad \Rightarrow \quad A = (\epsilon l) c$$

$$y = mx + c$$

$$\epsilon = \frac{m}{l}$$

$$\therefore c = \frac{A}{\epsilon l} \quad \Rightarrow \quad c = \frac{0.64}{1.366 \times 1}$$

$$= 0.344 \text{ M}$$




EXPERIMENT-5


AIM - To determine the strength of KMnO_4 soln with the help of potassium dichromate soln. using Mohr's salt as an intermediate.

MATERIAL REQUIRED - Mohr's salt, Potassium dichromate, n-phenyl antranilic acid, dil. sulphuric acid, KMnO_4 , Buoyette, Funnel, Stand, pipette, white tile, measuring flask, conical flask, test tubes.

THEORY - Acidic $\text{K}_2\text{Cr}_2\text{O}_7$ is a strong O.A. When it is added to ferrous ammonium sulphate soln containing dil. H_2SO_4 , only FeSO_4 is oxidised.



KMnO_4 being a powerful O.A. helps in estimating R.A. like Mohr's salt in acidic medium, according to -



PROCEDURE - 1) 10 ml of Mohr's salt soln was taken in a conical flask.

2) Half test tube of dil. H_2SO_4 was added to it.

3) The buvette was filled with known $\text{K}_2\text{Cr}_2\text{O}_7$ soln.

4) The colour of filtrate turned to green. On further addn. of $\text{K}_2\text{Cr}_2\text{O}_7$, it changed to

greenish purple.

- 5) 10 ml of Mohr's salt was taken in a conical flask and half test tube of dil. H_2SO_4 was added to it.
- 6) 3 concordant readings were taken.

RESULT - Strength of given $KMnO_4$ soln. is
1.007 g/L.

- PRECAUTION - 1) $KMnO_4$ is dark coloured, so read upper miniscus.
- 2) Don't use rubber cork burette.
 - 3) Use dil. H_2SO_4 for acidifying $KMnO_4$.

OBSERVATION - Table for known soln. -

SNO	Voln of soln taken	Burette Reading (ml)			concordant Reading (ml) Used
		Initial	Final	Total	
1.	10 ml	0	12.5	12.5	12.5 ml Diphenoxy Anisallic
2.	10 ml	13	25.5	12.5	12.5 ml
3.	10	26	38.5	12.5	12.5 ml

Table for unknown soln -

SNO	Voln of soln taken	Burette Reading			concordant Reading	Indicator Used
		Initial	Final	Total		
1.	10 ml	0	9.2	9.2	9.2 ml	KMnO ₄
2.	10 ml	10	19.2	9.2	9.2 ml	
3.	10 ml	20	29.2	9.2	9.2 ml	

CALCULATION -

$$N_1 V_1 = N_2 V_2 \quad \text{or, } N_1 = \frac{N_2}{40} \quad ; \quad V_1 = 12.5 ; V_2 = 10 \text{ ml}$$

$$\Rightarrow 10 \times N_2 = \frac{1}{40} \times 12.5$$

$$N_2 = 0.03 \text{ N}$$

$$N_1 V_1 = N_2 V_2$$

$$\Rightarrow N_1 \times 9.2 = 0.03 \times 10$$

$$N_1 = 0.03 \text{ N}$$

$$\text{Strength} = 0.03 \times \frac{158}{5} = 1.07 \text{ g/L}$$